## MERIT'S BLEEDBACK ALTERNATIVE

### INSTRUCTIONS FOR USE:

**DESCRIPTION:** MBA Hemostasis Valve with Integral Introducer PATENT PENDING

**INDICATIONS AND USAGE:** The MBA Hemostasis Valve is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a guiding catheter.

**CAUTION:** Federal (U.S.A.) Law restricts this device to use by or on the order of physician. The device is intended for single patient use only. Read instructions prior to use. Read manufacturer's instructions for the use of catheters, guide wires and introducers.

**INSTRUCTIONS FOR USE:** Inspect the device prior to use to verify that no damage has occurred during shipping.

1. Connect the sideport of the fully open hemostasis valve to the manifold assembly. Flush and fill the assembly with saline. To fill the valve section, ensure the thumb wheel is in the fully open position, open the slit valve by pushing in and slightly turning the integral introducer through the slit valve to fully open it. Place one finger over the male luer fitting and continue to fill the assembly until completely purged of air. Pull back to disengage the integral introducer and engage the secondary seal.

2. Connect the rotating male luer adapter of the hemostasis valve to the guiding catheter. Ensure the connection is airtight.

3. To open the valve for catheter or guide wire insertion, with the hemostasis valve in the fully open position, push in and slightly turn the integral introducer. The valve is now ready to accept any devices.

**Caution:** Valve will be fully open and allow for blood loss when integral introducer is engaged.

4. Insert the catheter or wire into the integral introducer and advance it to an appropriate distance into the vasculature. Disengage the integral introducer.

5. A fluid tight seal, providing little or no blood loss, will form around the interventional device yet does not restrict movement of the guide wire, catheter or introducer.

6. To secure the interventional device in place, turn the adjustable thumb wheel to tighten.

**Caution:** Do not inject any fluid if air bubbles are visible within the valve.

The narrowest rigid portion of the hemostasis valve has an inner diameter of 0.118 inches (0.305 cm) (approximately 9 Fr).

Manufacturer: Merit Medical Systems, Inc. South Jordan, Utah 84095  
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748 www.merit.com

Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland www.merit.com

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