



***Dual pressure valve:**
Standard pressure all the
time, over-pressurize when
needed with the click of a button.*

**PRESSURE
INFUSOR BAG**

Merit Medical Pressure Infusion Bag

CATALOG NUMBERS

PIB500 500ml Infusor Bag with
Dual Pressure Valve
12 per box, 48 per case

PIB1000 1000ml Infusor Bag with
Dual Pressure Valve
12 per box, 48 per case

MTM330E 3000ml Infusor Bag
2 per box, 10 per case

- *Two Pressure Settings*
- *Available in Multiple Sizes*
- *Loads Quickly and Easily*
- *Non-Sterile*



STANDARD PRESSURE
Valve is factory pre-set
for pressurization up to
300mmHg



OVER PRESSURE
Depress blue button to
set valve for over-pressure
mode allowing 300-
450mmHg for rapid
infusion



POLE HANGER
Easy and stable
placement onto IV pole



FLUID BAG HOOK
Bag attaches easily to rigid
hook which secures bag in
the event of pressure loss



ONE-HANDED PREP
Stopcock placed close
to bulb for easy inflation,
closure once fully pressur-
ized, and deflation



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

U.S. Patent No. 6,966,893

www.merit.com

400516005/A ID 040610