



Safety has its Advantages

*The next generation of temporary sharps holders starts
with the Merit Medical Shortstop Advantage.*

ShortStop[™] ADVANTAGE[™]

ShortStop^{ADVANTAGE}TM

TEMPORARY SHARPS HOLDERS

Merit Medical is pleased to introduce the ShortStop Advantage. This temporary sharps holder was designed to complement the classic ShortStop[®]. The ShortStop Advantage comes equipped with a new access point which reduces the possibility of accidental sharps injuries in the clinical setting.

ORDERING INFORMATION:

Catalog numbers:

SS10 – ShortStop[®]

SS20 – ShortStop AdvantageTM

Packaging: 25 per box - 4 boxes per case

FEATURES INCLUDE:


The new access point on the side of the device allows the user to place a fully loaded syringe with needle in a safe place avoiding any balance issues or bent needles.

A thick core resistant foam holds long needles and scalpels securely and minimizes particulate in needles.

Bright red container that meets OSHA standards and colorful targets for easy visualization.

An adhesive strip is attached to the bottom of the box to affix the ShortStop to the back table.

Thick ABS sides prevent sharps from protruding through the sides or bottom.



Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT • The Netherlands • +31 43 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115

www.merit.com

400777001/A ID 063006