


impress®

30CM ANGIOGRAPHIC CATHETERS


*Small additions to an
already impressive
catheter line*



impress®

30CM ANGIOGRAPHIC CATHETERS

Shorter catheter length enables the clinician to be closer to the working site on the forearm, while providing excellent torque.



OTHER PROCEDURAL SOLUTIONS



MERIT MAK™ & SMAK™ MINI ACCESS KITS

When little things mean a lot, get started with Merit.



PRELUDE® SHORT SHEATH INTRODUCERS

The shorter sheath provides access to dialysis fistulas and grafts. Can be used for declot, angioplasty, and thrombolysis procedures, as well as temporary hemodialysis.



SLIP-NOT® SUTURE RETENTION DEVICE

A unique and effective method to secure a purse-string suture following fistula intervention.



FOUNTAIN® INFUSION SYSTEMS

Exclusive gradient hole sizes allow uniform distribution of therapeutic agent along the entire length of the infusion segment.



PROGUIDE™ CHRONIC DIALYSIS CATHETER

A stepped tip dialysis catheter that's clearly a step up. Apposition bump deflection feature is designed to minimize recirculation and optimize flow rates.



SHORTSTOP® AND SHORTSTOP ADVANTAGE® TEMPORARY SHARPS HOLDERS

The ShortStop Advantage offers a side access feature which allows placement of a loaded syringe attached to a needle, thus preventing needle bending.

ORDERING INFORMATION

catalog number	french size	length	guide wire	tip shape	side holes
4303510PIG-NB	4F	30cm	.035" (0.89mm)	Pigtail	10
43035MODCB1	4F	30cm	.035" (0.89mm)	Cobra 1 mod	-
43035RIM	4F	30cm	.035" (0.89mm)	RIM	-
43035KA2	4F	30cm	.035" (0.89mm)	KA2	-
43035STR	4F	30cm	.035" (0.89mm)	Straight	-
43035BER	4F	30cm	.035" (0.89mm)	Berenstein	-
5303510PIG-NB	5F	30cm	.035" (0.89mm)	Pigtail	10
53035MODCB1	5F	30cm	.035" (0.89mm)	Cobra 1 mod	-
53035RIM	5F	30cm	.035" (0.89mm)	RIM	-
53035KA2	5F	30cm	.035" (0.89mm)	KA2	-
53035STR	5F	30cm	.035" (0.89mm)	Straight	-
53035BER	5F	30cm	.035" (0.89mm)	Berenstein	-
56535MRIMST	5F	65cm	.035" (0.89mm)	Mini RIM Short Tip	-

