



Renal Therapies & Accessories



CONTENTS

HEMODIALYSIS

- Prepare
- Access
- Navigate
- Implant
- Long-Term Access
- Treat

PERITONEAL DIALYSIS

- Prepare
- Place
- Select
- Integrate Assist
- Complete

Merit Medical's full suite of innovative dialysis products can help you improve your patient's quality of life when treating end-stage renal disease (ESRD). This comprehensive line of peritoneal dialysis and hemodialysis catheters, the HeRO[®] Graft, implantation kits, and accessories has been designed to meet the unique clinical needs of dialysis treatment.



HEMODIALYSIS

Prepare
Access
Navigate
Implant
Long-Term Access
Treat

Medallion®
Syringe




PAL™
Medication
Labeling System



Futura®
Safety Scalpel



FlowGuard™
Valved Introducer



HEMODIALYSIS
PREPARE

PROMOTE SAFE PREPARATION

Merit Medical has a complete line of products designed to promote the safety of clinicians and patients. **SHIELDED** scalpels, needles, and **SECURE** temporary sharps holders help protect clinicians from accidental sharps injuries. Multiple color options and smear-resistant labeling for syringes and containers help you prevent medication errors.

START WITH QUALITY


Achieve **QUICK** and effective access with our **SHARP** angiographic needles, scalpels, mini access kits designed to lessen trauma and reduce bleeding, and sheath introducers with excellent transitions for **SMOOTH** insertion.

HEMODIALYSIS
ACCESS

ShortStop Advantage®
Temporary Sharps Holder




Merit MAK™
Mini Access Kit




CentrosFLO®
Tunneler



ProGuide™
Tri-Ball Tunneler



Impress®
30 cm Diagnostic
Peripheral Catheter



Merit Laureate®
Hydrophilic Guide Wire



NAVIGATE WITH CONFIDENCE

LUBRICIOUS, hydrophilic guide wires reduce friction and enhance your ability to make rapid catheter exchanges. The **TAPERED** tips, kink resistant hub design, and unique braided shafts give Merit's diagnostic catheters both flexibility and **STEERABILITY**. Our diagnostic guide wires are manufactured with a unique **PTFE** precoating process which results in a consistently smooth surface for **SMOOTHER** tracking through difficult vasculature.

InQwire®
Diagnostic Guide Wire



CHOOSE PERFORMANCE

Merit's **STEP-TIP** catheter has a large internal lumen and apposition deflection feature to minimize recirculation and optimize flow. This catheter addresses the needs of clinicians using the **OVER-THE-WIRE** technique and those who prefer the valved, **PEEL AWAY SHEATH** insertion technique.

ProGuide™
Long-Term
Hemodialysis
Catheter

CentrosFLO®
Long-Term
Hemodialysis Catheter

FlowGuard™
Valved Introducer

FlowGuard™
Valved Introducer

IMPROVE PATENCY

CentrosFLO Long-Term Hemodialysis Catheters are the next generation, **SELF-CENTERING** split-tip dialysis catheters designed by Dr. Stephen Ash. The CentrosFLO's self-centering design was created to address the issue of fibrin sheath formation. The catheter has **SMOOTH** surfaces and two preshaped, **CURVED TIPS** to automatically "center" the catheter within the vein.

ProGuide™
Tri-Ball Tuner


CentrosFLO® Tuner

LONG-TERM ACCESS


ACCESS FOR HEMODIALYSIS

HeRO® (Hemodialysis Reliable Outflow) Graft is the **ONLY** fully **SUBCUTANEOUS** AV access solution clinically proven to maintain **LONG-TERM** access for hemodialysis patients with central venous stenosis, failing fistulas or catheter-dependency.

HeRO® Graft
Hemodialysis Reliable Outflow



Slip-Not®
Suture Retention
Device



basixTOUCH™
Inflation Device



EN Snare®
Endovascular
Snare System



ONE Snare®
Endovascular
Snare System



TREAT EFFECTIVELY

The Prelude Short Sheath provides single access to dialysis fistulas and grafts for **DECLOT**, angioplasty, thrombolysis and temporary hemodialysis. Two snare configurations give you the options you need to strip fibrin sheath formation. A complete line of therapeutic **INFUSION** systems allows you to achieve more uniform dispersion of therapeutic agents. Merit's latest innovation in inflation device technology offers increased volume and pressure capability for **INFLATION**. The Slip-Not **SUTURE RETENTION** device helps manage bleeding after an AV fistula intervention by securing purse-string sutures.


HEMODIALYSIS
TREAT

HEMODIALYSIS
TREAT


Squirt®
Fluid Dispensing
System



Fountain®
Infusion System




Prelude®
Short Sheath Introducer



PERITONEAL DIALYSIS

Prepare
Place
Select
Integrate
Assist
Complete



PREPARE

PLAN FOR OPTIONS

The **IMPLANTATION** stencil is designed to improve catheter placement **ACCURACY** by helping you select the optimal catheter style, size, implantation site and exit site for each patient. The stencil is conveniently included in every catheter kit and can be used for both left and right-sided placements. For upper abdomen and upper chest placements, Merit offers a complete set of three sterile **STENCILS** in each ExTended catheter kit.



Implantation Stencils

IMPLANT WITH SIMPLICITY

Merit offers a comprehensive implantation **KIT** for the single-site incision of a PD catheter, avoiding the need to place the catheter through an over-sized laparoscopic port. This complete set of components is designed to streamline the implantation procedure and promote **EFFICIENCY**. The kit includes a unique cuff implantor, dilators, trocar assembly, and a tunneler. Multiple kit **OPTIONS** are available to support your specific implantation approach and tunneling style.



Y-TEC® VP-211
Implantation System

PLACE

CHOOSE WITH CONFIDENCE

The Merit adult PD catheter has an internal diameter of 3.5 mm, allowing 35% **HIGHER FLOW** rates than other catheters currently on the market. Catheters are available in various **SHAPES**, and tip-to-cuff **LENGTHS** to accommodate most patient types including infant, pediatric, adolescent, and small, medium, and large adults. Our Classic and ARC™ shapes are designed to accommodate various tunneling techniques and exit site options. Our complete ExxTended™ catheter kit includes all of the necessary tools for upper abdomen and upper chest (pre-sternal) procedures.




Flex-Neck Classic
Peritoneal Dialysis Catheter

Flex-Neck ExxTended™
Peritoneal Dialysis Catheter



Flex-Neck® ARC™
Peritoneal Dialysis
Catheter



CUSTOMIZE FOR EFFICIENCY

Our **COMPLETE** ExxTended catheter kit includes all of the necessary tools for upper abdomen and upper chest (pre-sternal) procedures. The kit includes a measuring tape, a tunneling tool, three stencils for precise patient marking, and a measuring rod to gauge the rectus-to-skin distance for more **ACCURATE** catheter sizing and cutting.



FIT WITH PRECISION

Merit's PD accessories are specifically designed to **SUPPORT** catheter implantation, connection, and repair. The **REPAIR** Kit is designed to extend the life of the catheter. The metal Faller Trocar creates a clean, tight exit site, utilizing an antegrade tunneling approach. Three lengths of reusable, stainless steel implantation stylettes are specifically sized to **FIT** the interior diameter of our PD catheters.

Faller Trocar
Stainless Steel

Flex-Neck
PD Repair Kit


Catheter
Connectors

FINALIZE WITH CERTAINTY

The Embedding Tool is used for dialysis **PREPLANNING** in cases where residual kidney function exists or dialysis will be needed in the future. This tool allows the catheter to be embedded under the skin during the initial implantation procedure. The proximal end of the catheter can be explanted later for future use, once the patient is **READY** for dialysis. This allows the tissue ingrowth cuffs to adhere to the muscle and subcutaneous tissue **IN ADVANCE** of dialysis. Titanium connectors, plastic replacement connectors, and caps are available as needed by the PD clinic or when the catheter is placed.

Implantation
Stylette

Embedding® Tool



EFFICIENT. EFFECTIVE. SAFE.

AT EVERY STEP IN YOUR DIALYSIS
PROCEDURE, MERIT MEDICAL OFFERS
HIGH-QUALITY PRODUCTS DESIGNED
TO MEET YOUR CLINICAL NEEDS. FROM
PREPARATION TO INTERVENTION, OUR
PRODUCTS WILL HELP YOU ACHIEVE
THE CLINICAL, ECONOMICAL, AND
PROCEDURAL EFFICIENCY THAT IS
ESSENTIAL FOR DELIVERING THE
BEST POSSIBLE CARE.





Understand. Innovate. Deliver.™



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit.com

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg 8002
25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000

(Italiano)
+41 225180035

UK
0800 973 115