

MERIT **H₂**™ HYDROPHILIC GUIDE WIRE

SLIPPERY WHEN WET

GETS YOU THROUGH THE TOUGHEST CURVES



Slippery When Wet

HYDROPHILIC COATING

Increased lubricity enhances catheter tracking. Consistently smooth surface ensures integrity of the coating.

SOLID NITINOL CORE

Kink resistant core for precise manipulation and steerability in tortuous anatomy with 1:1 torque ratio.

INCREASED RADIOPACITY

Proprietary polymer technology provides enhanced visibility under fluoroscopy.

THROMBORESISTANT

A smooth polyurethane jacket and hydrophilic polymer significantly reduce the potential of thrombus formation on the guide wire.

STIFF SHAFT TECHNOLOGY

Provides strength and stability during catheter placement and exchange.

AVAILABILITY

Hydrophilic guide wires are available in a variety of diameters, lengths, and tip configurations to ensure the right wire is available for a wide range of diagnostic procedures.

Proprietary manufacturing process with strict process tolerances provide a dependable high performance wire.



ORDERING INFORMATION:

HYDROPHILIC STANDARD WIRE CATALOG NUMBERS

H2OSTDA3580	Standard	Angled	.035" (0.89mm)	80cm
H2OSTFA3580	Stiff	Angled	.035" (0.89mm)	80cm
H2OSTDS35150	Standard	Straight	.035" (0.89mm)	150cm
H2OSTDA35150	Standard	Angled	.035" (0.89mm)	150cm
H2OSTDA38150	Standard	Angled	.038" (0.97mm)	150cm
H2OSTFA35150	Stiff	Angled	.035" (0.89mm)	150cm

HYDROPHILIC SPECIALTY WIRE CATALOG NUMBERS

H2OSTDA35150LT	Standard	Angled	.035" (0.89mm)	150cm	Long Taper
H2OSTDA35180	Standard	Angled	.035" (0.89mm)	180cm	
H2OSTFS35180	Stiff	Straight	.035" (0.89mm)	180cm	
H2OSTFA35180	Stiff	Angled	.035" (0.89mm)	180cm	
H2OSTDA35260EX	Standard	Angled	.035" (0.89mm)	260cm	Exchange
H2OSTFS35260EX	Stiff	Straight	.035" (0.89mm)	260cm	Exchange
H2OSTFA35260EX	Stiff	Angled	.035" (0.89mm)	260cm	Exchange

HYDROPHILIC SMALL DIAMETER SPECIALTY WIRE CATALOG NUMBERS

H2OSTDA25150	Standard	Angled	.025" (0.64mm)	150cm	
H2OSTDA25180	Standard	Angled	.025" (0.64mm)	180cm	
H2OSTDS25180	Standard	Straight	.025" (0.64mm)	180cm	
H2OSTDA25260EX	Standard	Angled	.025" (0.64mm)	260cm	Exchange
H2OSTDA1880	Standard	Angled	.018" (0.46mm)	80cm	
H2OSTDA18150	Standard	Angled	.018" (0.46mm)	150cm	
H2OSTDA18180	Standard	Angled	.018" (0.46mm)	180cm	
H2OSTDA18260EX	Standard	Angled	.018" (0.46mm)	260cm	Exchange

All hydro wires are fixed core. Hydrophilic guide wires are sold 5 per box.

Torque devices may be purchased separately.

NOW NEW!
SMALL DIAMETER SPECIALTY WIRES



Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT-AIRPORT • The Netherlands • +31 - 43 - 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • **France** 0800 91 60 30 • **Germany** 0800 182 0871 • **Luxembourg** 0800 22522 • **Netherlands** 0800 022 8184 • **UK** 0800 973 115

www.merit.com

401814002/B ID 052505