




## KATETER UYUMLULUK TABLOSU

KATETER	İÇ ÇAP
Anjiyografik Kateter	0,046" (1168 µm)
Anjiyografik Kateter	0,040" (1016 µm)
EmboCath® Plus	0,028" (711 µm)
Merit Maestro®	0,027" (686 µm)
Merit Maestro®	0,024" (610 µm)
Merit Maestro®	0,020" (508 µm)

Uyumlu Partikül Büyüklüğü/Büyüklükleri



www.merit.com

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical EUROPE, THE MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • Hollanda • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, İrlanda • +353 (0) 91 703 733

Belirli ülkeler için ücretsiz telefon: Avusturya 0800 295 374 • Belçika 0800 72 906 (Flamanca) 0800 73 172 (Fransızca) • Danimarka 80 88 00 24 • Fransa 0800 91 60 30 • Finlandiya 0800 770 586

Almanya 0800 182 0871 • İrlanda (Cumhuriyet) 1800 553 163 • İtalya 800 897 005 • Lüksemburg 8002 25 22 • Hollanda 0800 022 81 84 • Norveç 800 11629 • İsveç 020 792 445 • Birleşik Krallık 0800 973 115

©2013 Merit Medical Systems, Inc. Tüm hakları saklıdır.

402840016/A ID 042813