


StayFIX[®]

FIXATION DEVICE



MERITMEDICAL[®] | Suture-free Fixation

StayFIX®

CLINICAL BENEFITS

Suture-Free - Noninvasive StayFIX means no painful sutures and no exposure to sharps or sharps-related injury.

Secure - Medical grade adhesive secures devices to patients for up to 7 days.

Easy - Patients can apply and remove the StayFIX by themselves.

Absorbent - Hydrocolloid material reduces skin irritation by absorbing fluid leaking from the insertion site.

Comfortable - Soft material conforms to the body and reduces catheter movement.

Water Resistant - Patients can shower without covering up the StayFIX or insertion site.

Other Merit drainage products



MAK-NV™
Introducer System



ReSolve® Locking
Drainage Catheter



ReSolve® Non-Locking
Drainage Catheter




ReSolve® Biliary Locking
Drainage Catheter



Merit Drainage Depot™
Drainage Bag

Application Instructions



Catalog Number	Fits Catheter Sizes	Box Quantity
680ME	5F -14F	25
685ME	12F-22F	20



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1-801-253-1600
1-800-35-MERIT
www.merit.com

Merit Medical Europe, Middle East,
and Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch)
0800 73 172 (French) | Denmark 80 88 00 24
France 0800 91 60 30 | Finland 0800 770 586 |
Germany 0800 182 0871 | Ireland (Republic) 1800 553 163
Italy 800 897 005 | Luxembourg 8002 25 22
Netherlands 0800 022 81 84 | Norway 800 11629
Sweden 020 792 445 | Switzerland 800 9000 5060
UK 0800 973 115

ID 080114 403035001/B